



iM20 Patient Monitor

Product Specifications

Product Specifications

NOTE:

The performance of the equipment with ☆ mark is determined to be essential performance.

A.1 Classification

| | |
|--------------------------|--|
| Anti-electroshock type | Class II equipment and internal powered equipment |
| Anti-electroshock degree | ECG (RESP), TEMP, IBP, C.O. CF SpO ₂ , NIBP, CO ₂ BF |
| Ingress Protection | IP44 (protected against splashing water and solid foreign objects ≥ 1.0 mm diameter) |
| Working system | Continuous operation equipment |
| Compliant with Standards | IEC 60601-1: 2005+A1 :2012; IEC 60601-1-2: 2014; EN 60601-1: 2006+A1 :2013; EN 60601-1-2: 2015; IEC 60601-2-49: 2011 |

A.2 Physical Specifications

| Product | Dimension | Max Weight | Comments |
|---------|-------------------------------------|------------|--|
| iM20 | 185 mm (W) ×116 mm (H) ×85.3 mm (D) | < 1.5 kg | Including battery, without accessories |
| EFM | 207 mm (W) ×116 mm (H) ×93.4 mm (D) | < 0.58 kg | Including sidestream CO ₂ module. |

A.3 Function Configuration

| Product | Standard Configuration | Optional Configuration |
|---------|--|--|
| iM20 | ECG (3- Electrode, 5- Electrode), RESP, SpO ₂ (EDAN), NIBP (EDAN), TEMP, Wi-Fi, USB interface | ECG (6- Electrode, 10- Electrode), SpO ₂ (Nellcor), NIBP (SunTech), IBP, C.O. |
| EFM | / | CO ₂ (Respironics LoFlo) |

Note: There is only one configuration for EFM module.

A.4 Environmental Specifications

The monitor may not meet the performance specifications given here if stored or used outside the specified temperature and humidity ranges.

When the monitor and related products have differing environmental specifications, the effective range for the combined products is that range which is common to the specifications for all products.

| | |
|-----------------------|------------------------------------|
| Main unit | |
| Temperature | |
| Working | +0 °C to +40 °C (32 °F ~ 104 °F) |
| Transport and Storage | -30 °C to +70 °C (-22 °F ~ 158 °F) |
| Humidity | |
| Working | 15%RH to 95%RH (non-condensing) |
| Transport and Storage | 15%RH to 95%RH (non-condensing) |
| Altitude | |
| Working | 61.5 kPa to 106 kPa |
| Transport and Storage | 61.5 kPa to 106 kPa |

| | |
|-----------------------|------------------------------------|
| EFM | |
| Temperature | |
| Working | +0 °C to +35 °C (32 °F ~ 95 °F) |
| Transport and Storage | -30 °C to +70 °C (-22 °F ~ 158 °F) |
| Humidity | |
| Working | 15%RH to 95%RH (non-condensing) |
| Transport and Storage | 15%RH to 95%RH (non-condensing) |
| Altitude | |
| Working | 61.5 kPa to 106 kPa |
| Transport and Storage | 61.5 kPa to 106 kPa |

NOTE:

The time required for the patient monitor to warm from the minimum storage temperature between uses until it is ready for intended use is at least 2 hours; the time required for the patient monitor to cool from the maximum storage temperature between uses until it is ready for intended use is at least 2 hours.

A.5 Out-Of-Hospital Transport Requirements

The monitor can be used in transport environments such as a road ambulance. For this purpose, the monitor meets the following requirements:

- EN 1789: 2007+A1: 2010 Road ambulances (Chapter 6 – Medical Devices).
- IEC/EN 60529 IP44 Specification for degrees of protection provided by enclosures.
- Radiated susceptibility 20 V/m according to Complies with ISO 80601-2-61: 2011. (SpO₂) and ISO 80601-2-55: 2011 (CO₂).

A.6 Power Supply

| | |
|---------|--------------------|
| Voltage | DC 11.1 V – 19.8 V |
| Current | 1.27 A – 2.3 A |

A.7 Display

| Display | Messages |
|---|--|
| Display screen: 5-inch color TFT, touch screen is configurable Resolution: 800 × 480 | One power on/off LED One battery charge LED One DC power LED One physiological alarm LED One technical alarm LED One alarm mute LED |

A.8 Battery

| | |
|----------------|---|
| Number | 1 |
| Battery Type | Lithium battery |
| Capacity | 11.1 V, 2400 mAh |
| Operating Time | 5.5 hrs (At 25 °C ± 2 °C, with (a) new fully charged battery, ECG (RESP)/TEMP/SpO ₂ module connected, NIBP automatic measurement mode at interval of 15 minutes, brightness set to “1”.) |
| Charging Time | ≤ 14 hrs (The monitor is on or in standby mode, 100% charge) |
| | ≤ 12.6 hrs (The monitor is on or in standby mode, 90% charge) |
| | 2.5 hrs (The monitor is off, 100% charge) |
| | 2.3 hrs (The monitor is off, 90% charge) |

A.9 Power Adapter

| | |
|------------------|--|
| AC Power Adapter | Input: 100 V - 240 Vac, 40VA; Output: 15 V \pm 5%dc, 24 VA |
| DC Power Adapter | Input: 12.4 V - 15.1 Vdc or 24.8 V - 30.3 Vdc, 1.6 A max; Output: 15 Vdc, 1 A max |

A.10 Data Management

Data Review

| | |
|-------------------------|---|
| Trend data | 1 hour, resolution: 1 s 150 hours, resolution: 1 min |
| Alarm events | Up to 200 sets |
| NIBP measurement data | 1200 sets |
| Arrhythmia events | Up to 200 sets |
| 12-Lead analysis result | Up to 50 sets |

Refer to Chapter *Review* for more information about data review.

Data Storage

A single piece of patient data maximally contains the following information:

| | |
|-----------------------------|--|
| Patient information | MRN, name, date of birth, date of admission, gender, type, height, weight, blood type, pace, doctor, bed No., department |
| Trend graph and trend table | 240 hours, resolution: 1 min |
| NIBP measurement review | 1200 sets |
| Alarm review | 200 sets |
| Arrhythmia event | 200 sets |
| 12-lead analysis review | 50 sets |
| Full disclosure waveforms | 3-lead/5-lead: 48 hours 12-lead: 35 hours |

Refer to Section *Storing Data in the Storage Device* for more information about storing data in the storage medium.

A.11 ECG

Complies with IEC 60601-2-25: 2011, IEC 60601-2-27: 2011.

| | |
|--|--|
| Lead Mode | <p>3 Electrodes: I, II, III</p> <p>5 Electrodes: I, II, III, aVR, aVL, aVF, V</p> <p>6 Electrodes: I, II, III, aVR, aVL, aVF, and leads responding to Va Vb.</p> <p>10 Electrodes: I, II, III, aVR, aVL, aVF, V1, V2, V3, V4, V5, V6</p> |
| Electrode Standard | AHA, IEC |
| ☆Display Sensitivity (Gain Selection) | 1.25 mm/mV ($\times 0.125$), 2.5 mm/mV ($\times 0.25$), 5 mm/mV ($\times 0.5$), 10 mm/mV ($\times 1$), 20 mm/mV ($\times 2$), 40 mm/mV ($\times 4$), AUTO gain |
| ☆Sweep | 6.25 mm/s, 12.5 mm/s, 25 mm/s |
| Bandwidth (-3dB) | <p>Diagnosis: 0.05 Hz to 150 Hz</p> <p>Diagnosis 1: 0.05 Hz to 40 Hz</p> <p>Monitor: 0.5 Hz to 40 Hz</p> <p>Surgery: 1 Hz to 20 Hz</p> <p>Enhanced: 2 Hz ~18 Hz</p> <p>Customized: High-pass Filter and Low-pass Filter (Refer to Section <i>Changing the ECG Filter Settings</i>)</p> |
| ☆CMRR (Common Mode Rejection Ratio) | <p>Diagnosis: > 95 dB</p> <p>Monitor: > 105 dB</p> <p>Surgery: > 105 dB</p> <p>Enhanced: > 105 dB</p> <p>Diagnosis 1: > 105 dB (when Notch is turned on)</p> <p>Customized: > 105 dB (Low-pass Filter < 40 Hz) > 95 dB (Low-pass Filter > 40 Hz)</p> |
| Hum Filter | In diagnosis, Diagnosis 1, monitor, surgery, enhanced modes: 50 Hz/60 Hz (Hum Filter can be turned on or off manually) |
| ☆ Differential Input Impedance | > 5 M Ω |
| ☆Input Signal Range | ± 10 mV PP |
| ☆Accuracy of Input Signal Reproduction | <p>An error of $\leq \pm 20\%$ of the nominal value of the output or ± 100 μV, whichever is greater.</p> <p>The total error and frequency response comply with IEC 60601-2-27: 2011, Sect. 201.12.1.101.1.</p> |

| | |
|--|--|
| ☆ Electrode Offset Potential Tolerance | ±500 mV |
| Auxiliary Current (Leads off detection) | Active electrode: < 100 nA Reference electrode: < 900 nA |
| ☆ Recovery Time After Defibrillation | < 5 s (measured without electrodes as IEC60601-2-27:2011, Sect. 201.8.5.5.1 requires.) |
| Leakage current of patient | < 10 µA |
| Scale signal | 1 mVPP, accuracy is ±5% |
| ☆System Noise | < 30 µVPP |
| ☆Multichannel Crosstalk | ≤ 5% of the input signal Complied with IEC 60601-2-27: 2011, Sect. 201.12.1.101.5. |
| ☆Frequency and Impulse Response | Frequency response: Input a 5 Hz, 1 mV sine wave signal, and the output signal amplitude remains within the range of 71% to 110% at 0.67 Hz and 40 Hz. Input a 1 Hz, 1.5 mV 200 ms triangular wave input signal, and the output shall be within 11.25 mm~15 mm. Impulse response: Displacement value: ≤ 0.1 mV Slope: ≤ 0.3 mV/s following the end of the pulse. Complied with IEC 60601-2-27: 2011, Sect. 201.12.1.101.8. |
| Sampling frequency | 1000 Hz |
| Sampling channel switch time | < 80 µS |
| A/D precision | 24 Bits (Minimum resolution: 0.07533 µV/LSB) |
| ☆ESU Protection | Cut mode: 300 W Coagulation mode: 100 W Restore time: ≤ 10 s |
| Electrosurgical Interference Suppression | Test according to ANSI/AAMI EC13:2002, Sect. 5.2.9.14. Complied with ANSI/AAMI EC13:2002, Sect. 4.2.9.14. |

| | |
|--|---|
| Pace Pulse | |
| Pulse indicator | <p>Pulse is marked if the requirements of IEC 60601-2-27: 2011, Sect. 201.12.1.101.12 are met:</p> <p>Amplitude: ± 2 mV to ± 700 mV</p> <p>Width: 0.1 ms to 2.0 ms</p> <p>Ascending time: 10 μs to 100 μs</p> |
| Pulse Rejection | <p>Pulse is rejected if the requirements of IEC 60601-2-27: 2011, Sect. 201.12.1.101.13 are met:</p> <p>Amplitude: ± 2 mV to ± 700 mV</p> <p>Width: 0.1 ms to 2.0 ms</p> <p>Ascending time: 10 μs to 100 μs</p> |
| Pace pulse detecting lead: one among I, II, III, aVR, aVL, aVF, V1 to V6 | |
| Minimum input slew rate (lead II) | > 2.5 V/S |
| ☆Baseline Reset Time | < 3 s |
| Heart Rate | |
| HR Calculation | |
| ☆Range | <p>ADU: 15 bpm to 300 bpm</p> <p>PED/NEO: 15 bpm to 350 bpm</p> |
| ☆Accuracy | $\pm 1\%$ or 1 bpm, whichever is greater |
| Resolution | 1 bpm |
| Sensitivity | ≥ 300 μ VPP |
| ☆QRS Detection Range | <p>The detection range has exceeded the requirement described in the standard:</p> <p>Width: 70 ms~120 ms for adult, 40 ms~120 ms for Pediatric/neonate.</p> <p>Amplitude: 0.5 mv~5 mv</p> <p>In adult mode, these two signals are not responded:</p> <ol style="list-style-type: none"> 1. when QRS amplitude of 0.15 mV or less is applied; 2. when QRS duration of 10 ms and QRS amplitude of 1 mV or less is applied. <p>Complied with IEC 60601-2-27: 2011, Sect. 201.12.1.101.15.</p> |

| | |
|------------------------------|---|
| PVC | |
| Range | ADU: (0 to 300) PVCs/ min PED/NEO: (0 to 350) PVCs/ min |
| Resolution | 1 PVCs/min |
| Pause/min | |
| Range | ADU/PED/NEO: (0 to 30) pauses/min |
| Resolution | 1 pause/min |
| ST value | |
| Range | -2.0 mV to +2.0 mV |
| Accuracy | -0.8 mV to +0.8 mV: ± 0.02 mV or 10%, whichever is greater. Beyond this range: not specified. |
| Resolution | 0.01 mV |
| HR Averaging Method | |
| Method 1 | Heart rate is computed by excluding the minimum and maximum values from the 12 most recent RR intervals and averaging the residual 10 RR intervals. |
| Method 2 | If each of three consecutive RR intervals is greater than 1200ms, then the four most recent RR intervals are averaged to compute the HR. |
| Range of Sinus and SV Rhythm | |
| Tachy | Adult: RR interval for 5 consecutive QRS complex ≤ 0.5 s. Pediatric/neonatal: RR interval for 5 consecutive QRS complex ≤ 0.375 s. |
| Normal | Adult: 0.5 s < RR interval for 5 consecutive QRS complex < 1.5 s. Pediatric/neonatal: 0.375 s < RR interval for 5 consecutive QRS complex < 1 s. |
| Brady | Adult: RR interval for 5 consecutive QRS complex ≥ 1.5 s. Pediatric/neonatal: RR interval for 5 consecutive QRS complex ≥ 1 s. |
| Range of Ventricular Rhythm | |
| V-Tach | 5 consecutive ventricular beats and ventricular HR ≥ 100 bpm. |

| | |
|---|---|
| Vent Rhythm | <p>Basic: 5 consecutive ventricular beats, and $40 \text{ bpm} \leq \text{ventricular HR} < 100 \text{ bpm}$.</p> <p>Advanced: 5 consecutive ventricular beats, and $20 \text{ bpm} \leq \text{ventricular HR} < 40 \text{ bpm}$.</p> |
| Vent Brady | <p>Basic: 5 consecutive ventricular beats, and $\text{ventricular HR} < 40 \text{ bpm}$.</p> <p>Advanced: 5 consecutive ventricular beats, and $\text{ventricular HR} < 20 \text{ bpm}$.</p> |
| Maximum start-up alarm time for Tachycardia | |
| V-Tach 1 mV 206bpm | <p>Gain 1.0: 10 s</p> <p>Gain 0.5: 10 s</p> <p>Gain 2.0: 10 s</p> |
| V-Tach 2 mV 195bpm | <p>Gain 1.0: 10 s</p> <p>Gain 0.5: 10 s</p> <p>Gain 2.0: 10 s</p> |
| Response time of Heart Rate Meter to Change in HR | <p>HR range: 80 bpm to 120 bpm Range : Within 11 s</p> <p>HR range: 80 bpm ~ 40 bpm Range : Within 11 s</p> |
| ☆Tall T-wave Rejection | Complied with IEC 60601-2-27: 2011, Sect. 201.12.1.101.17 minimum recommended 1.2 mV T-Wave amplitude |
| Accuracy of Heart Rate Meter and Response to Irregular Rhythm | <p>Complied with IEC 60601-2-27: 2011, Sect. 201.7.9.2.9.101 b) 4), the HR value after 20 seconds of stabilization is displayed as follows:</p> <p>Ventricular bigeminy: $80 \text{ bpm} \pm 1 \text{ bpm}$</p> <p>Slow alternating ventricular bigeminy: $60 \text{ bpm} \pm 1 \text{ bpm}$</p> <p>Rapid alternating ventricular bigeminy: $120 \text{ bpm} \pm 1 \text{ bpm}$</p> <p>Bidirectional systoles: $91 \text{ bpm} \pm 1 \text{ bpm}$</p> |
| Time to Alarm for Heart Rate alarm conditions | <p>Asystole alarm: $\leq 10 \text{ s}$</p> <p>HR low alarm: $\leq 10 \text{ s}$</p> <p>HR high alarm: $\leq 10 \text{ s}$</p> |

| | | | |
|---|----------------------------------|----------------|-------------------|
| Arrhythmia analyses | Asystole | V-Fib/V-Tach | Couplet |
| | Vent Rhythm | PVC Bigeminy | PVC Trigeminy |
| | Tachy | R on T | PVC |
| | Irr Rhythm | Brady | Missed Beat |
| | Pacer not Pacing | Vent Brady | Pacer not Capture |
| | VEB | Run PVCs | Acc. Vent Rhythm |
| | IPVC | Non-Sustain VT | Multiform PVCs |
| | Pauses/min High | Pause | Afib |
| | PAC Bigeminy | PVCs High | Low Voltage(Limb) |
| | ExtremeBrady | PAC Trigeminy | Wide QRS Tachy |
| | Sustain VT | ExtremeTachy | V-Tach |
| 12-lead ECG Synchronization Analysis | Average parameters of heart beat | | |
| | Heart rate (bpm) | | |
| | Time limit of P wave (ms) | | |
| | PR interval (ms) | | |
| | QRS interval (ms) | | |
| | QT/QTc (ms) | | |
| | P-QRS-T AXIS | | |

A.12 RESP

| | |
|---------------------------------|--|
| Method | Impedance between RA-LL, RA-LA |
| Measurement lead | Options are lead I and II. The default is lead II. |
| Calculation Type | Manual, Automatic |
| Respiration excitation waveform | Sinusoid, 45.6 kHz ($\pm 10\%$), $< 500 \mu\text{A}$ |
| Measuring Sensitivity | Within baseline impedance range: 0.3Ω |
| Waveform bandwidth | 0.2 Hz to 2.5 Hz (-3 dB) |
| Baseline Impedance Range | 200 Ω to 2500 Ω (leads cables 1 K Ω resistance) |
| ☆RR Measuring Range | |

| | |
|--------------------------|---|
| ☆ Adult | 0 rpm to 120 rpm |
| ☆ Neo/Ped | 0 rpm to 150 rpm |
| Resolution | 1 rpm |
| ☆ Accuracy | |
| ☆ Adult | 6 rpm to 120 rpm: ± 2 rpm 0 rpm to 5 rpm: not specified |
| ☆ Neo/Ped | 6 rpm to 150 rpm: ± 2 rpm 0 rpm to 5 rpm: not specified |
| ☆ Gain Selection | $\times 0.25$, $\times 0.5$, $\times 1$, $\times 2$, $\times 3$, $\times 4$, $\times 5$ |
| ☆ Sweep | 6.25 mm/s, 12.5 mm/s, 25 mm/s |
| ☆ Apnea Alarm Time Setup | 10 s, 15 s, 20 s, 25 s, 30 s, 35 s, 40 s; default value is 20 s. |

A.13 NIBP

Complies with IEC 80601-2-30: 2009+A1: 2013.

EDAN Module

| | |
|---|--|
| Technique | Oscillometry |
| Mode | Manual, Auto, Continuous, Sequence |
| Measuring interval in AUTO Mode(unit: minute) | 1/2/3/4/5/10/15/30/60/90/120/180/240/360/480 and User Define (default is 2.5) |
| Continuous | 5 min, interval is 5 s |
| Measuring Parameter | SYS, DIA, MAP, PR |
| Pressure Unit | kPa, mmHg, cmH ₂ O |
| ☆ Measuring Range (Applicable for CE registration area) | |
| ☆ Adult Mode | SYS: 25 mmHg to 290 mmHg DIA: 10 mmHg to 250 mmHg MAP: 15 mmHg to 260 mmHg |
| ☆ Pediatric Mode | SYS: 25 mmHg to 240 mmHg DIA: 10 mmHg to 200 mmHg MAP: 15 mmHg to 215 mmHg |
| ☆ Neonatal Mode | SYS: 25 mmHg to 140 mmHg DIA: 10 mmHg to 115 mmHg MAP: 15 mmHg to 125 mmHg |

| | |
|---|--|
| ☆Measuring Range (Applicable for FDA registration area) | |
| ☆Adult Mode | SYS: 40 mmHg ~ 270 mmHg DIA: 10 mmHg ~ 215 mmHg MAP: 20 mmHg ~ 235 mmHg |
| ☆Pediatric Mode | SYS: 40 mmHg ~ 230 mmHg DIA: 10 mmHg ~ 180 mmHg MAP: 20 mmHg ~ 195 mmHg |
| ☆Neonatal Mode | SYS: 40 mmHg ~ 135 mmHg DIA: 10 mmHg ~ 100 mmHg MAP: 20 mmHg ~ 110 mmHg |
| ☆Alarm Type | SYS, DIA, MAP |
| ☆Cuff Pressure Measuring Range | 0 mmHg ~ 300 mmHg |
| Pressure Resolution | 1 mmHg |
| ☆Maximum Mean Error | ±5 mmHg |
| ☆Maximum Standard Deviation | 8 mmHg |
| Maximum measuring period | |
| Adult/Pediatric | 120 s |
| Neonate | 90 s |
| Typical Measuring Period | 20 s to 35 s (depend on HR/motion disturbance) |
| Dual Independent Channel Overpressure Protection | |
| Adult | (297±3) mmHg |
| Pediatric | (245±3) mmHg |
| Neonatal | (147±3) mmHg |
| Pre-inflation Pressure | |
| Adult Mode | 80 mmHg, 100 mmHg, 120 mmHg, 140 mmHg, 150 mmHg, 160 mmHg, 180 mmHg, 200 mmHg, 220 mmHg, 240 mmHg Default: 160 mmHg |
| Pediatric Mode | 80 mmHg, 100 mmHg, 120 mmHg, 140 mmHg, 150 mmHg, 160 mmHg, 180 mmHg, 200 mmHg Default: 140 mmHg |
| Neonatal Mode | 60 mmHg, 70 mmHg, 80 mmHg, 100 mmHg, 120 mmHg Default: 100 mmHg |

| | |
|-----------------------|---|
| Venipuncture pressure | |
| Adult | Default: 60 mmHg Options: 20 mmHg, 30 mmHg, 40 mmHg, 50 mmHg, 60 mmHg, 70 mmHg, 80 mmHg, 90 mmHg, 100 mmHg, 110 mmHg, 120 mmHg |
| Pediatric | Default: 40 mmHg Options: 20 mmHg, 30 mmHg, 40 mmHg, 50 mmHg, 60 mmHg, 70 mmHg, 80 mmHg |
| Neonatal | Default: 30 mmHg Options: 20 mmHg, 30 mmHg, 40 mmHg, 50 mmHg |

SunTech Module

| | |
|--|---|
| Method | Oscillometric |
| Mode | Manual, Auto, Continuous and Sequence |
| Measuring Interval in AUTO Mode (unit: minute) | 1/2/3/4/5/10/15/30/60/90/120/240 and User Define |
| ☆Measuring Parameter | SYS, DIA, MAP, PR |
| ☆Measuring Range | |
| ☆Adult Mode | SYS: 40 mmHg ~ 260 mmHg DIA: 20 mmHg ~ 200 mmHg MAP: 26 mmHg ~ 220 mmHg |
| ☆Pediatric Mode | SYS: 40 mmHg ~ 230 mmHg DIA: 20 mmHg ~ 160 mmHg MAP: 26 mmHg ~ 183 mmHg |
| ☆Neonatal Mode | SYS: 40 mmHg ~ 130 mmHg DIA: 20 mmHg ~ 100 mmHg MAP: 26 mmHg ~ 110 mmHg |
| ☆Alarm Type | SYS, DIA, MAP |
| Pressure Resolution | 1 mmHg |
| ☆Maximum mean error | ±5 mmHg |
| ☆Maximum standard deviation | 8 mmHg |
| Maximum measuring period | |
| Adult | 130 s |

| | |
|-------------------------|---|
| Adult (Sports Mode) | 120 s |
| Pediatric | 90 s |
| Neonate | 75 s |
| Overpressure protection | |
| Adult/Pediatric | < 300 mmHg |
| Neonate | < 150 mmHg |
| Pre-inflation Pressure | |
| Adult Mode | 120 mmHg, 140 mmHg, 150 mmHg, 160 mmHg, 180 mmHg, 200 mmHg, 220 mmHg, 240 mmHg, 260 mmHg, 280 mmHg Default: 160 mmHg |
| Pediatric Mode | 80 mmHg, 100 mmHg, 120 mmHg, 140 mmHg, 150 mmHg, 160 mmHg, 180 mmHg, 200 mmHg, 220 mmHg, 250 mmHg Default: 140 mmHg |
| Neonatal Mode | 60 mmHg, 70 mmHg, 80 mmHg, 90 mmHg, 100 mmHg, 120 mmHg, 140 mmHg Default: 90 mmHg |

A.14 SpO₂

Complies with ISO 80601-2-61: 2011.

EDAN Module

| | |
|---------------------|---|
| Measuring Range | 0% to 100% |
| Alarm Range | 20% to 100% |
| Resolution | 1% |
| ☆Data Update Period | 1 s |
| ☆Accuracy | |
| ☆Adult /Pediatric | ±2% (70% to 100% SpO ₂) |
| | Undefined (0% to 69% SpO ₂) |
| ☆Neonate | ±3% (70% to 100% SpO ₂) |
| | Undefined (0% to 69% SpO ₂) |
| Sensor | |
| Red light | (660±3) nm |
| Infrared light | (905±10) nm |

| | |
|----------------------|------------------------------|
| Emitted light energy | < 15 mW |
| PI | |
| Measuring Range | 0-10, invalid PI value is 0. |
| Resolution | 1 |

Nellcor Module

| | | |
|----------------------|---|--------------------------------------|
| Measuring Range | 1% to 100% | |
| Alarm Range | 20% to 100% | |
| Resolution | 1% | |
| ☆ Data Update Period | 1 s | |
| ☆ Accuracy | DS-100A, OXI-A/N(Adult) D-YS (Adult and Pediatric) OXI-P/I (Pediatric) | ± 3% (70% to 100% SpO ₂) |
| | MAX-A, MAX-AL, MAX-N, MAX-P, MAX-I, MAX-FAST (Adult and Pediatric) | ±2% (70% ~ 100% SpO ₂) |
| | MAX-A, MAX-AL, MAX-N, MAX-P, MAX-I, MAX-FAST (Adult and Pediatric) | ±3% (60% ~ 80% SpO ₂) |
| | If sensor is used for neonate as recommended, the accuracy will be larger greater than adult by ±1. | |
| Sensor | Wave length: approximately 660 nm and 900 nm | |
| | Emitted light energy: < 15 mW | |

NOTE:

Information about the wave length range can be especially useful to clinicians (for instance, when photodynamic therapy is performed).

A.15 PR

| | | Measuring range | Accuracy | Resolution |
|------------------------|---------|-------------------|-----------------------------|------------|
| PR (SpO ₂) | EDAN | 25 bpm to 300 bpm | ±2 bpm | 1 bpm |
| | Nellcor | 20 bpm to 300 bpm | ± 3 bpm (20 bpm to 250 bpm) | 1 bpm |

| | | Measuring range | Accuracy | Resolution |
|-----------|---------|-------------------|--|------------|
| PR (NIBP) | EDAN | 40 bpm to 240 bpm | ± 3 bpm or 3.5%, whichever is greater | 1 bpm |
| | SunTech | 30 bpm to 220 bpm | ± 3 bpm or $\pm 2\%$, whichever is greater | 1 bpm |
| PR (IBP) | EDAN | 20 bpm to 300 bpm | 30 bpm to 300 bpm: ± 2 bpm or $\pm 2\%$, whichever is greater; 20 bpm to 29 bpm: undefined | 1 bpm |

A.16 TEMP

Complies with ISO 80601-2-56: 2009.

| | |
|-------------------------|---|
| Channel | 2 |
| Sensor type | YSI-10K and YSI-2.252K |
| Technique | Thermal resistance |
| Measuring Mode | Direct Mode |
| Position | Skin, oral cavity, rectum |
| Measuring Range | 0 °C to 50 °C(32 °F to 122 °F) |
| Resolution | 0.1 °C (0.1 °F) |
| ☆Accuracy ¹ | ± 0.3 °C |
| Refresh Time | Every 1 s to 2 s |
| Temperature Calibration | At an interval of 5 minutes to 10 minutes |
| Transient Response Time | ≤ 30 s |

Note 1: The accuracy consists of two parts, as following:

- Accuracy (not including sensor): ± 0.1 °C
- Sensor accuracy: $\leq \pm 0.2$ °C

A.17 IBP

Complies with IEC 60601-2-34: 2011.

| | |
|------------------|-----------------------------|
| Technique | Direct invasive measurement |
| ☆Measuring Range | |
| Art | 0 mmHg to + 300 mmHg |

| | | |
|-----------------------------------|---|---------------------|
| PA | -6 mmHg to + 120 mmHg | |
| CVP/RAP/LAP/ICP | -10 mmHg to + 40 mmHg | |
| P1/P2 | -50 mmHg to + 300 mmHg | |
| Resolution | 1 mmHg | |
| ☆ Accuracy (not including sensor) | ± 2% or ±1 mmHg, whichever is greater ICP: 0 mmHg to 40 mmHg: ±2% or ±1 mmHg, whichever is greater; -10 mmHg to -1 mmHg: undefined | |
| Pressure Unit | kPa, mmHg, cmH ₂ O | |
| Pressure sensor | | |
| Sensitivity | 5 (μV/V/mmHg) | |
| Impedance | (300 to 3000) Ω | |
| Filter | DC~ 12.5 Hz; DC~ 40 Hz | |
| Zero | Range: ±200 mmHg | |
| Pressure Calibration Range | IBP (excluding ICP) | 80 mmHg to 300 mmHg |
| | ICP | 10 mmHg to 40 mmHg |
| Volume displacement of MSI | 7.4 x 10 ⁴ mm ³ /100mmHg | |

A.18 CO₂

Complies with ISO 80601-2-55: 2011.

| | | |
|-------------------------|--|--------|
| Applicable Patient Type | Adult, pediatric and neonatal patients | |
| Technique | Infra-red Absorption Technique | |
| Unit | mmHg, %, Kpa | |
| ☆ Measuring Range | | |
| ☆ EtCO ₂ | 0 mmHg ~ 150 mmHg | |
| ☆ FiCO ₂ | 3 mmHg ~ 50 mmHg | |
| ☆ AwRR | 2 rpm ~ 150 rpm (Sidestream) | |
| Resolution | EtCO ₂ | 1 mmHg |
| | FiCO ₂ | 1 mmHg |
| | AwRR | 1 rpm |

| | |
|--------------------------------------|--|
| ☆EtCO ₂ Accuracy | ±2 mmHg, (0 to 40) mmHg |
| | ±5% of reading, (41 to 70) mmHg |
| | ±8% of reading, (71 to 100) mmHg |
| | ±10% of reading, (101 to 150) mmHg |
| | ±12% of reading, RR is over 80 rpm (sidestream) |
| ☆AwRR Accuracy | ±1 rpm |
| Operation Mode | Measure, standby |
| Sample Gas Flowrate (sidestream) | (50 ±10) ml/min |
| O ₂ Compensation | |
| Range | 0%~ 100% |
| Resolution | 1% |
| Default | 16% |
| Barometric pressure compensation | User setup |
| Anesthetic Gas Compensation | |
| Range | 0% ~ 20% |
| Resolution | 0.1% |
| Default | 0.0% |
| Balance Gas Compensation | Room air, N ₂ O, helium |
| Stability | |
| Short Term Drift | Drift over 4 hours < 0.8 mmHg |
| Long Term Drift | 120 hours |
| Total System Response Time | 4.7 s |
| Alarm Type | EtCO ₂ , FiCO ₂ , AwRR |
| Apnea Alarm Delay | 10 s, 15 s, 20 s, 25 s, 30 s, 35 s, 40 s; default value is 20 s. |
| Data Sample Rate | 100 Hz |
| Sensor Response Time (sidestream) | < 3 seconds, including transport time and rise time |

Interfering Gas and Vapor Affects on EtCO₂ Measurement Values:

| Gas or vapor | Gas level (%) | Quantitative effect/Comments |
|--------------|---------------|------------------------------|
|--------------|---------------|------------------------------|

| | | |
|---------------|----|--|
| Nitrous oxide | 60 | Dry and Saturated Gas |
| Halothane | 4 | (0 ~ 40) mmHg: ± 1 mmHg additional error |
| Enflurane | 5 | (41 ~ 70) mmHg: $\pm 2.5\%$ additional error |
| Isoflurane | 5 | (71 ~ 100) mmHg: $\pm 4\%$ additional error |
| Sevoflurane | 5 | (101 ~ 150) mmHg: $\pm 5\%$ additional error |
| Xenon | 80 | *Additional worst case error when compensation for P_B , O_2 , N_2O , anesthetic agents, or helium is correctly selected for the actual fractional gas constituents present. |
| Helium | 50 | |
| Desflurane | 15 | Desflurane: The presence of desflurane in the exhaled breath at concentrations greater than 5% will positively bias Carbon Dioxide values by up to an additional 3 mmHg at 38 mmHg. Xenon: The presence of Xenon in the exhaled breath will negatively bias Carbon Dioxide values by up to an additional 5 mmHg at 38 mmHg. |

Barometric Pressure on EtCO₂ Measurement Values:

| |
|--|
| Quantitative effect |
| Ambient Barometric, Operational |
| (0 ~ 40) mmHg: ± 1 mmHg additional error |
| (41 ~ 70) mmHg: $\pm 2.5\%$ additional error |
| (71 ~ 100) mmHg: $\pm 4\%$ additional error |
| (101 ~ 150) mmHg: $\pm 5\%$ additional error |
| *Additional worst case error when compensation for P_B , O_2 , N_2O , anesthetic agents, or helium is correctly selected for the actual fractional gas constituents present. |

NOTE:

Respiration Rate accuracy was verified by using a solenoid test setup to deliver a square wave of known CO₂ concentration to the device. 5% and 10% CO₂ concentrations were used. Respiration rate was varied over the range of the device. Pass/Fail criteria was comparison of the respiratory rate output from the sensor to the frequency of the square wave.

A.19 C.O.

| | |
|--------------------|--------------------------|
| Technique | Thermodilution Technique |
| Measure Parameters | C.O., TB, TI |
| Measuring Range | |

| | |
|------------|---|
| C.O. | 0.1 L/min to 20 L/min |
| TB | 23 °C to 43 °C (73.4 °F to 109.4 °F) |
| TI | -1 °C to 27 °C (30.2 °F to 80.6 °F) |
| Resolution | |
| C.O. | 0.1 L/min |
| TB, TI | 0.1 °C (+0.1 °F) |
| Accuracy | |
| C.O. | ± 5% or ± 0.2 L/min, whichever is greater |
| TB | ± 0.1 °C (not including sensor) |
| TI | ± 0.1 °C (not including sensor) |

NOTE:

At least 90% of the C.O. data should reside inside the bounded region, and the lower 95% confidence interval should not exceed 85%.

A.20 Wi-Fi

| | |
|---------------------------------------|--|
| IEEE | 802.11b/g/n |
| Frequency | 2.4 GHz ISM band |
| Modulation | OFDM with BPSK, QPSK, 16-QAM, and 64-QAM 802.11b with CCK and DSSS |
| Typical Transmit Power (± 2 dBm) | 17 dBm for 802.11b DSSS 17 dBm for 802.11b CCK 15 dBm for 802.11g/n OFDM |

A.21 Interface

A.21.1 USB Interface

| | |
|--------------------------|-------------------------------|
| Number of USB Interfaces | Standard: 1 |
| Drive Mode | OTG, USB1.0/2.0 protocol |
| Power Supply | 5 VDC $\pm 5\%$, 150 mA Max. |
| Interface Type | Micro USB-type port |