

SWING-ARM ELITE CARE CLINER

6 SERIES

6940 (STANDARD) / 6950 (X-LARGE)



Ease Patient Transfers Without Sacrificing Patient Comfort...

- Dual swing-arms with foot-activated release levers
- Dual side-tables
- Built-in headrest with cover
- Trendelenburg positioning
- 5" Total-lock and directional-lock casters

ENHANCING THE PATIENT & CAREGIVER EXPERIENCE



SPECIFICATIONS

SWING-ARM ELITE CARE CLINER	6940	6950
DESCRIPTIONS	STANDARD	X-LARGE
Seat Height	21" (53.3 cm)	
Seat Width	21" (53.3 cm)	25" (63.5 cm)
Seat Depth	20" (51 cm)	
Overall Height	48.5" (123.2 cm)	
Overall Width*	33.5" (85.1 cm)	37.5" (95.3 cm)
Range of Seat Angle	8° - 32°	
Range of Back Angle	62° - 38°	
Back Height above Seat	32.5" (82.6 cm)	
Overall Length (Back Upright)	41" (111.8 cm)	
Overall Length (Back Reclined)	76.5" (194.3 cm)	
Weight	104 lbs (47.2 kg)	127 lbs (57.6 kg)
Weight Capacity	350 lbs (158.7 kg)	450 lbs (204 kg)

* Add 9" to overall chair width for each raised side-table.


SUGGESTED OPTIONS & ACCESSORIES



HEAT
(H)



HEAT & MASSAGE
(H/M)



IV POLE
Left, Right, or
Dual

OVER 50
DESIGNER
COLORS
AVAILABLE

Additional options and accessories available.

Please visit our website to view a full selection of options or contact your local representative.

WARRANTY Contact Winco or visit our website for warranty details.

COLORS Visit our website or contact Winco for available color selections and/or swatches

NOTE Winco reserves the right to change specifications without notice.

CERTIFICATIONS 1) Winco recliners meet ANSI/BIFMA X5.4 for reliability testing 2) Winco foam meets California Technical Bulletin 117, Flammability Test Procedure for Resilient Filling Materials used in Upholstered Furniture 3) NR or TB options are available. NR meets AB2998 regulations and TB meets NFPA 101 and ASTM E1537-16 Standards for Upholstered Furniture 4) Heat and Massage, Heat Only Options: Meets electrical safety standards, UL 1647 and CAN/CSA-C22.2 No.68-92



CONTACT US

To arrange a demonstration or to speak with one of our associates, please feel free to contact us at:

800-237-3377 | WincoMfg.com | 5516 SW 1st Lane, Ocala, FL 34474 | Info@WincoMfg.com