

# **Advantages**

### **EMR** Connectivity

Scale has the ability to interface with a computer, EMR software or other electronic device via various Health o meter® Professional Connectivity Solutions

Scale functions include the ability to calculate the patient's Body Mass Index (BMI)

## Low Profile Platform

Low platform height requires only minimal step up, offering patient a more comfortable and stable weighing experience

### **▼ Value-Added Features**

Operating functions (LB/KG Conversion, LB/KG Lock Out, BMI, Zero Out/Tare, Hold/Release, Auto Zero, Auto Off) included at no additional cost

# **Compact Size**

The slim design of the scale and waisthigh display require a limited amount of office space



# $\boldsymbol{\omega}$ S S



Model	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
499KL	26.625" x 17.375" x 11.75" / 3.15	22 lb	892076002527







• Capacity: 500 lb / 220 kg

• Resolution: 0.2 lb / 0.1 kg

Platform Size: 13." (w) x 16." (d) x 23/8" (h) / 349 mm (w) x 419 mm (d) x 60 mm (h)

Height Rod: Optional 499KLROD

• Display: 1" LCD

• Wheels: 2

Connectivity: USB

 Power Source: 6 AA batteries included, Optional 120V adapter (order # ADPT31)

 Functions: LB / KG Conversion, LB / KG Lock Out, Body Mass Index (BMI), Zero Out / Tare, Hold / Release, Auto Zero, Auto Off

• 2 Year Limited Warranty



Health o meter® Professional Scales 9500 West 55th Street USA McCook, IL 60525-7110 Phone: 1-800-815-6615

Fax: 1-800-654-7330